

Find out what is in your food – FAST!

Yolanda Fintschenko, Ph.D.
Formerly Thermo Fisher Scientific,
Director of New Technology and Business Development,
LabSmith
Livermore, CA

Why Speed Matters



FOR IMMEDIATE RELEASE – Sun Prairie, WI - January 22, 2010 – The **Wisconsin Cheeseman**[®] announced today that it is **recalling cheese log/cheese ball** products in conjunction with the **voluntary recall** initiated by **Parkers Farm, Inc.** on January 5, 2010. Parkers Farm, Inc., located in Coon Rapids, Minnesota, believes some of its food items have the potential to be contaminated with **Listeria monocytogenes bacteria**. The recall was a result of a sampling done by the state of Wisconsin and the state of Minnesota



Speed Matters for Residue Analysis



U.S. Food and Drug Administration


FOR IMMEDIATE RELEASE
June 28, 2007

FDA News

FDA Detains Imports of Farm-Raised Chinese Seafood Products Have Repeatedly Contained Potentially Harmful Residues

During targeted sampling...FDA repeatedly found that farm – raised seafood imported from China were contaminated with antimicrobial agents that are not approved for this use in the United States.

Automate What? The Analysis Workflow



Sample Prep Automation

- Requires a good understanding of sample variability
- Focuses on removing human error and saving time
- May require an improved sample prep workflow in order to automate

Automation is NOT always robotics

Automating QuEChERS?


- 15 steps HOMOGENIZE SAMPLE
- 2.5 hours for 12 – 40 samples
- Skilled Technician

- Consumables EXTRACTION
 - Tubes
 - Salts and solvents


- Other equipment CLEANUP
 - Vortexer
 - Centrifuge

SOLVENT EXCHANGE

Automated Solid Phase Microextraction (SPME)



SPME with Quantum GC



Benefits of automation and manual sample preparation


SPME	Manual sample prep
Extraction and concentration in one step	Extraction and concentration in different steps- more time
Minimum use of solvent (or none)	Use of more solvent
Headspace and liquid extraction	Liquid extraction
Less maintenance for source & column	More matrix contamination
Works well for low fat matrices	Works for high and low fat samples
Excellent selectivity & sensitivity w/ TSQ Quantum GC	

In-Line Sample Prep Automation – Turbulent Flow Chromatography

- Turbulent flow chromatography

Size

- Turbulent flow chromatography
- Velocity dependent
- Molecular weight dependent



Cost/Benefit Analysis TurboFlow Technology

Costs Include

Initial investment for TurboFlow system

Service maintenance

TurboFlow column

Benefits Include

Simplify sample preparation by eliminating SPE/PPT/LLE

Direct injection of fluid samples

Reduce ion suppression

Can be multiplexed to increase sample throughput

Cost saving in time, labor and solvent

Solving Salmonella in Peanut Butter



Peanut Butter Salmonella Outbreak Rages On

CDC: Infected People Ate Salmonella-Contaminated Peanut Butter Crackers

By [Daniel J. DeNoon](#)

WebMD Health News

Jan. 21, 2009 -- More than 125 consumer peanut butter products, from more than 70 companies, have been recalled in the ongoing U.S. [salmonella](#) outbreak.

The most recent person to get sick fell ill on Jan 8...

Challenges for Analysis of Pathogens

- Multiple serotypes
- Traceability
- Enrichment
- Complying with regulations
- Identification
- Time




Salmonella typhimurium



Salmonella enteritidis plate


The Sample Prep Bottleneck



The Modular Approach – Functional Fittings




Fluidic and electronic breadboard



Zero dead volume valves and connectors




**Luer Lock
Capillary connectors**




Refillable functional cartridges

Automating Centrifugation



- Chirica, *et al.* Anal. Chem. **2006**, **78**, 5362-5368

Room for R&D: Automated Pathogen Detection



The automated sample preparation (ASP) system had a total cycle times of 15 minutes for *Bacillus subtilis*.

Vandernoot *et al*, Anal. Chem. 2007, 79, 5763-5770

Could this apply to assays for food?

FOOD
navigator.com

■ Portable labs detect water toxins faster

By Wai Lang Chu, 29-Sep-2005

...designed and built at Sandia/California. *"This on-site monitoring approach would enhance current utility monitoring systems that require water samples to be sent to laboratories for analysis, which sometimes takes days for results,"*

Conclusions

- Automation is more than robotics
- Advantages include
 - Precision
 - People free to solve complex problems
 - Reduced Time
 - Reduced Cost
- There is room for automation in food safety analysis!

Acknowledgements

- Thermo Fisher
 - Dr. Stuart Cram
 - Dr. Guiping Lu
 - Dr. James Chang
 - Dr. Kimihide Ohmichi
 - Dr. Michal Godula
 - Vinny Paez
 - Gerry Broski
- LabSmith
 - Monica Gottero
 - Kirsten Pace
- Sandia National Laboratories
 - Dr. Gabriela Chirica
 - Dr. Victoria Vandernoot
 - Ron Renzi
 - Dr. Julie Fruetel
 - Brent Haroldsen




Thank You!

www.labsmith.com

yfintschenko@labsmith.com

Seven Classes of Vet Drugs in Honey

- Macrolides (2)
- Amphenicols (2)
- β -Lactams (8)
- Aminoglycosides (3)
- Acaricides (12)
- Sulfonamides (17)
- Tetracyclines (5)



ES Antimicrobial
ES Antimicrobial

Common Approaches to Rapid Pathogen Analysis in Food

- Immunoassay
 - Test strips
 - Latex agglutination assays/reverse passive latex agglutination assays
 - ELISA/ELFA
 - Immunomagnetic
 - Immunodiffusion
- Nucleic acid
 - PCR
 - RT-PCR
- Biochemical
- ATP