


T116-A360 AND T116-A132 ADAPTER DIMENSIONS

INCHES
(MM)




MATERIAL: PEEK

PART NUMBER	COMPONENT SIZE [MALE]	TUBING SIZE [FEMALE]
T116-A360	1/16"	360 μ m
T116-A132	1/16"	1/32"



SEE C360-100 / A132-100 DRAWING FOR DETAILED DIMENSIONS OF THIS COMPONENT





SECTION B-B
SCALE 4 : 1



T116-A132 ADAPTER HAS SCRIBE LINES TO DISTINGUISH IT FROM T116-A360



SECTION A-A
SCALE 4 : 1



6-32UNC