

360 μ m CAPILLARY AND 1/32" INTERCONNECT DIMENSIONS

INCHES
(MM)

MATERIAL: ULTEM OR OPTICALLY CLEAR ACRYLIC

DRAWING SHOWN IS OF C360-204 CROSS. DIMENSIONS CALLED OUT ARE THE SAME FOR ALL C360 AND T132 INTERCONNECTS. T132 INTERCONNECTS ARE MARKED 1/32" AT EACH PORT




TABLE 1	
DIM -A-	MATERIAL
.150" [3.81]	ULTEM
.17" [4.3]	ACRYLIC

