


# LUER-LOCK ADAPTOR DIMENSIONS

**INCHES  
(MM)**

MATERIAL: PEEK



T132-300 ADAPTER HAS  
SCRIBE LINE TO DISTINGUISH  
FROM OTHER SIZES

TABLE 1

ADAPTOR P/N	DIM -A-	DIM -B-	ADAPTOR SIZE
C360-300	0.550 [14.0]	2-56UNC	360µm CAPILLARY
T132-300	.550 [14.0]	2-56UNC	1/32" TUBING
T116-300	.626 [15.9]	6-32UNC	1/16" TUBING

DIM -B-  
SEE TABLE 1

